## **Steps For Using Physical Baxter**

## Boot Up

- 1. Ensure Ethernet is unplugged from Baxter
- 2. Reboot Laptop
- 3. Turn Baxter on (white push button on back)
- 4. Wait for following
  - a. Green LED ring lights up on Baxter's head
  - b. Ultrasonic sensors turn on (listen for them)
  - c. No more load screen on face screen
- 5. Plug Ethernet into Baxter
- 6. Run "cd baxter\_ws → ./bax.sh → source devel/setup.bash"
- 7. Run "rostopic list" and make sure all of baxter's topics show up
- 8. Good to run your code

## Power Down

- 1. End your code (ctrl C)
- 2. Run "rosrun baxter tools tuck arms.py -t"
- 3. Run "rosrun baxter tools enable robot.py -d"
- 4. Turn Baxter off (press white push button on back DO NOT HOLD)